www.tiger-x.com

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Our ceramic crucibles are made of first refractory material and do not contain any impurity which could pollute or alter the alloys.

- for a longer life of the ceramic crucibles we would suggest you to preheat them before their utilization, in the same furnace where the casting cylinders are preheated, for 5-10 minutes at least and to remove, after the melting process the slag from the crucible.
- Note that many alloys while melting cover themselves with a superficial oxidations layer and falsify completely the reading of the temperature causing an overheating of the alloy which will disintegrate the crucible.
This problem is easy to resolved by breaking the layer with a ceramic stirring rod.
- The smokes produced by the alloys may make dirty optic system of the reading device. In this case the alloy as well as the crucible will be overheated in a noxious way.
Always keep the optical reader clean.
- Introducing pieces of metal in such a size that they could jammed between the sides walls of the crucibles and could cause the shattering of the crucible. The pieces while heating, expand and practice a considerable pressure on the side-walls and cause shatter of the crucible.

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